# R2P SlenGuide™

## Guiding Catheter / Cathéter Guide / Catéter Guía / Cateter Guia

**CAUTION:** Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.  
**ATTENTION:** Conformément au droit fédéral des États-Unis, ce dispositif ne peut être vendu que par ou sur prescription d'un médecin.  
**ATENCIÓN:** La Ley Federal (EE. UU.) restringe la venta de este dispositivo a un médico o bien por orden de este.  
**ATENÇÃO:** A Lei Federal dos E.U.A. restringe a venda deste dispositivo por médicos ou por sua prescrição.

- Stacking limit by #  
Ne pas empiler plus de #  
Aplazamiento máximo de #  
Não empilhar a mais de # alturas
- Fragile  
Fragil  
Fragil
- Keep dry  
Conserver au sec  
Mantener seco  
Mantner seco
- Usable length  
Longueur  
Longitud  
Comprimento utilizável
- Keep away from sunlight  
Crainit la chaleur  
Mantener alejado de la luz solar  
Mantner ao resguardo da luz solar

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


TERUMO CORPORATION  
 44-1, 2-CHÔME, HATAGAYA, SHIBUYA-KU, TOKYO 151-0072, JAPAN MADE IN JAPAN

TERUMO CORPORATION 2019-08

## ENGLISH

**DESCRIPTION**  
 The R2P SlenGuide consists of a Guiding Catheter and an Inner Guide, which, when combined, are used as a dual catheter system. The distal part of the Guiding Catheter is coated with a hydrophilic polymer, which becomes highly lubricious when wet.

- Instruction for use  
**CAUTION:** Federal law (USA) restricts this device to sale by or on the order of a physician.  
**<Specifications of Guiding Catheters>**  
 • Usable length of Guiding Catheter: 150cm, 120cm  
 • Hydrophilic polymer coating: 300 mm  
 • Outer diameter of Guiding Catheter: 7Fr. (2.37 mm)  
 • Inner diameter of Guiding Catheter: 0.087" (2.20 mm)  
**<Specifications of Inner Guide>**  
 • Inner diameter of Inner Guide: 0.045" (1.15 mm)  
 • Inner Guide protruding length: 30 mm

Please read all instructions prior to use

**INDICATIONS FOR USE**  
 The R2P (radial to peripheral) SlenGuide is designed for the introduction of interventional and diagnostic devices for the peripheral vasculature of the lower extremities.


- CONTRAINDICATIONS**
- Patients with a contraindication to antiplatelet and/or anticoagulation therapy.
  - Pregnant patients or patients who may be pregnant.
  - Patients with known upper extremity vascular disease, extreme tortuosity, anomalous radial artery takeoff, severe subclavian stenosis or severe atherosclerosis.
  - Patients with known Buerger's disease or Raynaud's Phenomenon.
  - Presence of arterio-venous fistula.
  - Patients with known anatomical abnormalities of the upper extremity/artery.
  - Patients with known dissecting thoracoabdominal aortic aneurysm.
  - Do not use within the cerebral vasculature, the carotid arteries and the coronary arteries.
  - Do not use via femoral access site.

- Do not use any drugs containing organic solvents and/or oleaginous contrast media. Such agents may cause damage to the SlenGuide catheter.
- Do not use the Inner Guide alone. It should always be used together with the Guiding Catheter. The SlenGuide Inner Guide is not designed to be used alone.
- Use the Guiding Catheter together with the Inner Guide when inserting the Guiding Catheter or adjusting the position of the Guiding Catheter.
- Do not inject drugs or contrast media using power injection as this may cause rupture of the thin wall of the SlenGuide.
- Do not expose to or wipe this product with organic solvents or alcohol for disinfection. This may result in possible separation or breakage of the product or abrasion of the hydrophilic coating.
- The SlenGuide Guiding Catheter contains a metallic core, do not use with any inappropriate equipment (e.g. MRI).
- Prior to beginning radial artery access, conduct a screening test such as an Allen test to ensure the radial access is appropriate for the patient.
- The SlenGuide should only be used for transradial access by a physician who is well trained in introduction of interventional or diagnostic devices via transradial access.
- If introducing devices via a transradial approach, standard departmental radial artery protocol should be followed to prevent radial artery occlusion, injury or spasm.


**PRECAUTIONS FOR STORAGE**  
 Avoid exposure to water, direct sunlight, extreme temperature or humidity during storage.

**DIRECTIONS FOR USE**

- 1. Prepare to assemble**
- 1-1. Carefully remove the holder tubes containing the Guiding Catheter and Inner Guide from the package.  
**Precautions**
    - Do not bend this product at the edge of the holder tube, as this may result in possible separation or breakage of the product.
    - Do not use if the unit package or the product has been damaged or soiled.
  - 1.2. Flush the inside of Guiding Catheter by injecting heparinized physiological saline solution through the Guiding Catheter hub, using a syringe (see Fig. 4).


- 1-3. Also flush the inside of Inner Guide by injecting heparinized physiological saline solution through the Inner Guide hub, using a syringe (see Fig. 5).


- Precaution**  
 When gently inserting the Inner Guide into the Guiding Catheter, avoid trapping air bubbles into the combined device system.

- 2-3. Connect a hemostasis valve to the Inner Guide hub to reduce possible risk of blood backflow (see Fig.10).

- 3. Remove this product from the holder tube**
- 3-1. Gently remove the Inner Guide/Guiding Catheter combination from the holder tube.  
**Precaution**
    - If you feel any resistance, stop manipulating the Guiding Catheter. Flush it by injecting more heparinized physiological saline solution through the holder hub using a syringe and keep the surface of the Guiding Catheter wet when removing from the holder tube.






- 4. Insert this product into the blood vessel**
- 4-1. Use a cut down, Seldinger or percutaneous technique with an introducer sheath of an appropriate size to introduce the Inner Guide/Guiding Catheter combination into the vessel.
  - 4-2. Insert a guidewire of an appropriate size, previously immersed in heparinized physiological saline solution, into the introducer sheath.
  - 4-3. After that, insert the Inner Guide/Guiding Catheter combination into the introducer sheath over the guidewire in the vessel.  
**Precautions**
    - Carry out the appropriate anticoagulant and antiplatelet therapy during operation.
    - If you feel any resistance, stop manipulating the Inner Guide/Guiding Catheter combination.
    - Be careful not to damage the distal end of the guidewire when advancing the guidewire into the introducer sheath, as this may cause damage to the Inner Guide/Guiding Catheter combination to the introducer sheath and/or to the interventional/diagnostic devices (see Fig. 11).

- 5. Advance to the target position**
- 5-1. Advance the Inner Guide/Guiding Catheter combination through the introducer sheath over the guidewire into the target vessel under fluoroscopy.  
**Precautions**
    - Manipulate the Inner Guide/Guiding Catheter combination slowly and carefully not to damage the hydrophilic coating on the distal part of the Guiding Catheter.
    - When advancing the Inner Guide/Guiding Catheter combination into the vessel, be careful not to cause any damage to the vessel wall.
    - The Inner Guide/Guiding Catheter combination may occlude the artery. Care must, therefore, be taken to avoid complete obstruction of blood flow.

- Under fluoroscopy, be careful to lead the Guiding Catheter tip right into the desired position, to avoid any damage to or perforation of the vessel.

1-4. Flush the outside of the Guiding Catheter by injecting heparinized physiological saline solution through the holder hub, using a syringe (see Fig.6).

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TERUMO, R2P et SlenGuide are trademarks of TERUMO CORPORATION.