TERUMO


2020-06 Rev.5

Glidecath™ XP

Angiographic Catheter
Cathéter d' Angiographie
Angiographiekatheter
Catéter Angiográfico

Catetere Angiografico
Angiografie Catheter
Angiografisk kateter
Cateter Angiográfico

RC34E103-03

REF
Catalogue number
Référence du catalogue
Bestellnummer
Número de catálogo
Codice prodotto
Catalogo number
Artikelnummer
Referência

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3. Gain access to the artery by using a percutaneous or cut down technique.
CAUTION Consider the use of systemic heparinization.
4. Insert the guide wire alone into the artery. Then proceed to advanced the catheter into the artery over the guide wire.
CAUTION To avoid damage to the catheter after it has been advanced into the vessel, manipulate the guide wire carefully, particularly when negotiating a bend in the catheter and/or when passing through the catheter's tip.
5. Manipulate the catheter slowly and carefully in the artery.
WARNING Never advance or withdraw an intraluminal device against resistance until the cause of resistance is determined by fluoroscopy. Failure to exercise proper caution may result in damage to the vessel or catheter. Separation of the catheter may occur requiring retrieval in some cases.
6. When the catheter tip has reached the branch of the desired vessel, remove the guide wire through the catheter.
7. While confirming the location of the catheter tip under fluoroscopy, advance the catheter to the desired site and perform angiography.
WARNING Before starting infusion, verify that the catheter has not been kinked or blocked. Failure to abide by this warning may cause the catheter to break/rupture/separate, resulting in damage to the vessel.
CAUTION When injecting contrast media, do not exceed the maximum injection pressure of 1,000 psi (6.895 kPa).
8. When using the catheter to guide a micro catheter, insert the micro catheter carefully into the lumen of the catheter.
CAUTION Use a micro catheter that is below the maximum guide wire size compatible to this catheter.
9. After the procedure is completed, draw the catheter back from the site. Insert the guide wire into the catheter until it extends slightly beyond the distal end of the catheter. Carefully remove the catheter and guide wire together.

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ENGLISH

INSTRUCTION FOR USE
The GLIDECATH XP should be used by a physician who is well trained in manipulation and observation under fluoroscopy.
Rx ONLY CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.
Read all instructions prior to use. To avoid complications, observe all warnings and precautions throughout these instructions.

PACKAGING / STORAGE
This product has been sterilized by ethylene oxide gas. Sterile and non-pyrogenic in an unopened and undamaged packaging.
CAUTIONS
• Do not use if the unit package or the product has been damaged or soiled. Use immediately after opening the unit package. Dispose of safely after single use to avoid risk of infection.
• For single use only. Do not reuse. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
• Do not store at extreme temperatures and humidity. Avoid direct sunlight.

INDICATION FOR USE
The GLIDECATH XP is intended for use in angiographic procedures. It delivers radiopaque media and therapeutic agents to selected sites in the vascular system. It is also used to lead a guide wire or a catheter into the target site.

DESCRIPTION
Angiographic catheter with hydrophilic coating.
WARNINGS AND PRECAUTIONS
WARNINGS Failure to observe these warnings may result in damage to the vessel, abrasion of the hydrophilic coating, damage to or breakage/separation of the catheter that may necessitate retrieval.
• Do not heat or bend the catheter tip. Abrasion of the hydrophilic coating or damage to the catheter may result.
• Consider the use of systemic heparinization to prevent or reduce the possibility of thrombus formation on the surface of the catheter.
• The surface of the distal portion of the catheter must be completely wet with heparinized saline solution to remain lubricated. Proper caution must be taken to wet the catheter prior to and during use.
• Never advance the guide wire briskly and/or force it into the catheter when the catheter is bent or twisted. This may cause breakage/separation of the catheter, resulting in damage to the vessel.


PRECAUTIONS
• Select the catheter of optimal tip shape and size, taking into account the site in which it is to be advanced to, as well as the patient's anatomy.
• When using a drug or a device with the GLIDECATH XP, the operator should have a complete understanding of the properties/ characteristics of the drug or device and exercise due caution to avoid damage to the catheter.
• Do not exceed the maximum injection pressure. Contrast media should be injected at 37°C.

CONTRAINDICATIONS
In certain medical conditions, angiography may be contraindicated or, special precautionary measures may need to be taken prior to and/or during the procedure. The following is a list of some of these medical conditions:
• The acute phase of myocardial infarction • A serious arrhythmia • Serious serum electrolyte imbalance • An allergy to or reaction to the contrast medium • Renal dysfunction • Blood coagulopathies • Congested Heart Failure • Some Respiratory Disorders
• Mental Disease • Pregnancy

COMPLICATIONS
Catheterization may be accompanied by but not limited to the following:
• Artery dissection • Artery injury • False aneurysm • Artery perforation • Arteriovenous fistula • Spasm • Vascular thrombosis
• Distal embolism • Hematoma • Bradycardia • Hemorrhage • Infection and pain at the puncture site

DIRECTIONS FOR USE
1. Carefully open the sterile pouch and gently remove the catheter from the package. Flush it by injecting heparinized saline solution through the catheter hub using a syringe. Soak the catheter prior to use in heparinized saline to lubricate the surface.
CAUTIONS
• Do not use if the catheter has been damaged or any other anomaly is observed.
• Prior to use make sure the catheter is immersed in heparinized saline solution. Keep the surface of the catheter wet during use.
• Some catalogue numbers have a metal shape keeper located in the distal portion of the catheter. Remove it before wetting the surface with heparinized saline solution.
• Employ an aseptic technique during removal from the package and use.

2. Insert a guide wire of appropriate size into the catheter through its hub and advance the wire to approximately 5cm beyond the catheter's distal tip.
CAUTION The GLIDECATH XP can accommodate a maximum guide wire size of 0.038" (0.97 mm) diameter.

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©TERUMO CORPORATION 2020-06
All brand names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners.