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English / French / Spanish / Portuguese

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
NaviCross®

Support Catheter
Cathéter de soutien
Catéter de apoyo
Cateter de Apoio

Icons and text for REF (Catalogue number), LOT (Batch code), # (Contents), STERILE EO (Sterilized using ethylene oxide), and various warnings like 'Keep dry', 'Fragile', 'Keep away from sunlight', 'P max', 'Non-pyrogenic'.

TERUMO CORPORATION
44-1, 2-CHOME, HATAGAYA, SHIBUYA-KU, TOKYO 151-0072, JAPAN
MADE IN JAPAN
©: Registered Trademark
©: Marque Enregistrée
©: Marca Registrada
©: Marca registrada

TERUMO

© Terumo Corporation 2019-04

ENGLISH

Please read all instructions prior to use.

INDICATIONS

NaviCross 0.018" is intended to guide and support a guide wire during access of the peripheral vasculature, allow for wire exchanges and provide a conduit for the delivery of saline or diagnostic contrast agents.

CONTRAINDICATIONS
Patients with a contraindication to anti-platelet and/or anti-coagulation therapy.
Women who are pregnant or are supposed to be pregnant.
Do not use within the cerebral vasculature and the coronary arteries.

COMPLICATIONS
Possible complications associated to catheterization include, but are not limited to, the following:
Bleeding complications - Allergies to drugs - Distal embolization - Hypotension
Infection and pain at puncture site - Arteriovenous fistula - Hematoma
Bradycardia and Palpitation - Intra-vessel thrombosis
Femoral pseudoaneurysm/Pseudoaneurysm formation - Arterial spasm
Arterial injury/perforation/dissection - Arterial embolism/occlusion - Arrhythmia
Hemorrhage and hemorrhagic shock - Fever/chill - Nausea and vomiting

WARNINGS AND PRECAUTIONS

WARNINGS
Do not heat or bend the catheter tip. Abrasion of the hydrophilic coating or damage to the catheter may result.
Carefully handle the product under fluoroscopy. If any resistance is felt while handling the product, immediately stop the manipulation and find out the cause to avoid damage to blood vessels and separation or breakage of the product.
The entire procedure should be carried out aseptically.
Before inserting/withdrawing the product, clean the surface of the guide wire with gauze moistened with saline solution. Advancing/withdrawing the product over a guide wire with residual blood on its surface or a guide wire which is not fully wet may result in separation or breakage of the product.
Do not torque the product excessively while the distal part of the product crosses the stenosis or is in the stent.
Careful when inserting/withdrawing the product through an opening of the stent struts to avoid damage to the product. A scratch by the stent struts may result in separation or breakage of the product.
Careful when exchanging a guide wire while leaving the product in the vessel. Carefully insert a guide wire into the product. If any resistance is felt, stop the manipulation and remove the product together with the guide wire to avoid separation or break of the product.
Perform appropriate anticoagulant or antiplatelet therapy according to the patient's condition to avoid complications, such as thrombotic embolization.

PRECAUTIONS
This product must be used by physicians who are well trained in manipulation. When using a drug or a device with the product, the operator should have a complete understanding of the properties/characteristics of the drug or device and exercise due caution to avoid damage to the catheter.
Consider the use of systemic heparinization.
Infusion flow rates for this product at maximum injection pressures for saline and iohexol (10.6 mPa a 37°C) are provided below. Do not exceed the maximum permissible injection pressure. Contrast media should be injected at 37°C.

Table with 3 columns: Guide Wire compatibility / Shaft diameter, Flow Rate (Saline, Iohexol), Maximum Permissible Pressure (300 psi / 2,068 kPa).

This product has been sterilized by ethylene oxide gas. For single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
Do not torque the product excessively if it is already bent, to avoid separation or breakage of the product.
Sterile and non-pyrogenic in an unopened and undamaged unit packaging. Do not use if the unit package or the product has been damaged or soiled.
Use immediately after opening the package. Dispose of safely after single use to avoid risk of infection.
Use the product prior to the expiry date indicated on the package.
The catheter should be examined to verify functionality and ensure that its size and shape are suitable for the specific procedure for which it is to be used.

DIRECTIONS FOR USE

1. Carefully remove the product in its holder from the package and carefully take the product out of the holder. Flush it by injecting heparinized saline solution through the catheter hub using a syringe. Soak the catheter prior to use in heparinized saline to lubricate the surface.

CAUTIONS
Do not bend the product at the edge of the holder. The product may break or separate.
Do not use if the catheter has been damaged or any other anomaly is observed.
Prior to use make sure the catheter is immersed in heparinized saline solution. Keep the surface of the catheter wet during use.
Apply an aseptic technique during removal from the package and use.

2. Insert the product over the wire into the hemostasis valve attached to the guiding catheter or sheath.

CAUTIONS
Take care not to damage the product, when a guide wire is inserted from its distal end.
Use a guide wire of an appropriate diameter.
Use a guiding catheter or sheath with an appropriate inner diameter.
Make sure that the hemostasis valve is open enough for insertion of the product. If not, the product may break or separate.
Do not advance the guide wire briskly and/or force it into the catheter when the catheter is bent or twisted. This may cause breakage/separation of the catheter, resulting in damage to the vessel.

3. Advance the guide wire and the product to the target vessel under fluoroscopy.

CAUTIONS
Manipulate the product slowly and carefully in the vessel because of its hydrophilic coating.
Do not insert the product into a vessel with a smaller diameter than the outer diameter of the product.
Do not torque the product if it is or seems stuck in order to avoid separation or breakage of the product.

4. In case injection of contrast media is required, withdraw the guide wire and inject the contrast media from the catheter hub.

CAUTIONS
Before starting infusion, verify that the catheter has not been kinked or blocked. Ignoring this warning may cause the catheter to break/rupture/separate, resulting in damage to the vessel.
When injecting contrast media, do not exceed the maximum permissible pressure of 300psi (2,068kPa).

5. Open the hemostasis valve. Withdraw the product along the guide wire leaving the guide wire inside the vessel. Close the hemostasis valve after removing the product.

CAUTIONS
Confirm position of the guide wire under fluoroscopy when the product is removed.
Remove the product, the guide wire and the guiding catheter or sheath together if any resistance is felt while withdrawing the product.
After removal, the catheter should be rinsed with heparinized saline solution to remove blood residues from the catheter's surface.
If blood residues do not come off, wipe gently once with a gauze soaked with heparinized saline solution. Do not use disinfectants, which can compromise the integrity of the hydrophilic coating. If necessary, flush the lumen of the catheter with heparinized saline solution to remove blood residues.

PRECAUTIONS FOR STORAGE

Avoid exposure to water, direct sunlight, extreme temperature or high humidity storage.

CAUTION

Federal law (USA) restricts this device to sale by or on the order of a physician.

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