TERUMO
Nagare™ Bi-Directional Steerable Sheath Instructions for Use
Device Description
The Nagare™ Bi-Directional Steerable Sheath consists of a dilator and steerable sheath which are designed to provide intracardiac access and flexible, stable catheter positioning in the cardiac anatomy. A hemostasis valve in the handle minimizes blood loss during catheter introduction and/or exchange. A side port with 3-way stopcock is provided for air or blood aspiration, fluid infusion, blood sampling and pressure monitoring. A deflection knob on the handle deflects the tip of the sheath clockwise ≥ 180° and counterclockwise ≥ 90°. Distal tip vent holes facilitate aspiration and a radiopaque marker provides visualization of the sheath tip fluoroscopically.

Indications
The Nagare™ Bi-Directional Steerable Sheath is indicated for introducing various cardiovascular catheters into the vasculature and into the chambers of the heart including the left side of the heart through the interatrial septum.

Contraindications

- Previous intra-atrial septal patch
- Known or suspected atrial myxoma
- Myocardial infarctions within the last two weeks
- Unstable angina
- Recent Cerebral Vascular Accident (CVA)
- Patients who do not tolerate anticoagulation therapy
- Patients with a known or suspected allergy to contrast media
- Patients with an active infection
- Presence of atrial thrombus

Warnings

- Do not alter this device in any way
- Only those physicians who are trained in transseptal procedures and delivery systems should use this device
- Do not reuse this device. Thorough cleaning of biological and foreign material is not possible. Adverse patient reactions may result from reuse of this device.
- Maintain continuous hemodynamic monitoring throughout the procedure
- Always observe acceptable hemodynamics prior to advancing the dilator or any other component
- Always monitor device movement during insertion. Do not advance or withdraw against resistance as this may result in patient injury or death, or damage to the device.
- Always withdraw components slowly to minimize the vacuum created during withdrawal
- Always aspirate slowly to minimize the vacuum created during aspiration
- Aspirate all air via the side port prior to fluid infusion
- Provide continuous heparinized saline infusion while the sheath remains in the vessel
- Fibrin may accumulate in or on the sheath tip during the procedure. To prevent dislodgment of potential thrombus, aspirate when removing dilator or catheter.
- Prior to removing the steerable sheath, reinsert the guidewire through the sheath, reintroduce the dilator over the guidewire, straighten the steerable sheath, then remove the dilator, guidewire, and sheath as a unit
- Maximum in-vivo time: 7 hours

Precautions

- Carefully reading the Instructions for Use of this device and any associated transseptal device will help reduce the potential risks and complications associated with the transseptal technique such as air emboli and/or perforation of the aorta and left atrium or the vasculature.
- Inspect all components before use. Do not use if the package or items in the kit appear to be damaged or defective.
- The French size specified represents the inner diameter of the sheath.
- Do not attempt to insert a catheter having a distal tip or body size larger than the sheath size indicated on the package label
- Do not attempt to use a guidewire larger than the maximum diameter specified on the package label
- Prior to inserting the device into the patient, pre-assemble the steerable sheath and dilator
- During handling and insertion, use caution not to create excessive bends in this device
- Frequently aspirate and saline flush the sheath to minimize the potential for thrombus formation
- Do not remove dilator or catheter rapidly
- Catheters should be inserted and removed through the sheath when the sheath is deflected to 90° or less
- If resistance is met when advancing or withdrawing guidewire or sheath, determine the cause and correct before continuing with the procedure
- Aspirate slowly, only from the side port
- Inject or saline flush only from the side port
- Discard all used sheaths, catheters and sterile accessories in accordance with hospital procedures
- Certain conditions may require special consideration when using this product. These may be, but are not limited to Enlarged Aortic Root, Small Left Atrium, Marked Distortion of the Thorax (Configuration (i.e. Kyphosis or Scoliosis).
- Store at ambient conditions.
- This device is intended for one time use only. Do not reprocess or reuse. Reuse can cause patient injury and/or the communication of infectious disease(s) from one patient to another, or device malfunction.

Potential Complications

The following complications may occur during the use of this device, but are not limited to:

- | | | |
|---------------------|------------------------------|------------------------|
| • Air embolism | • Thrombus formation | • Kidney damage |
| • Cardiac tamponade | • Allergic reaction | • Blood loss |
| • Infection | • Arrhythmia | • Atrial septal defect |
| • Vascular damage | • Vasovagal reaction | • Death |
| • Hematoma | • Ischemia | |
| • Perforation | • Stroke or systemic embolus | |

Please consult the adverse events associated with the use of the cardiovascular catheters.

How Supplied

The Nagare™ Bi-Directional Steerable Sheath is supplied sterile and packaged within a single pouch. The package contents are listed on the respective pouch and box labels.

Directions for Use

Note: the Nagare™ Bi-Directional Steerable Sheath is designed for temporary (<7 hours) use within the body.

Carefully read the Instructions for Use of this device to reduce potential risks associated with transseptal technique such as air emboli or perforation of the aorta or left atrium. Only advance transseptal needle when stylette is in place. Only physicians specially trained in transseptal procedures should use this device in a procedure involving access through the interatrial septum. Fluoroscopy should be used to confirm positioning throughout the procedure. Transseptal procedures should be performed only in facilities appropriately equipped and staffed to perform such procedures. Lab capabilities include but are not limited to:

- Systemic pressure monitoring
- Contrast media injection and management of reactions to contrast media
- Intracardiac pressure monitoring equipment
- Pericardiocentesis
- Surgical backup
- Anticoagulation therapy and monitoring
- Vital sign monitoring throughout the procedure

Device Preparation

- Prior to removing the sheath and accessories from packaging, inspect to verify the integrity of the seal. Use sterile techniques when handling this device.
- Prior to inserting the device into the patient, test the deflection mechanism to ensure it is operational. To deflect the tip of the sheath, rotate the deflection knob in a clockwise or counterclockwise direction. Wet the dilator shaft with sterile saline and insert the dilator into the hemostasis valve of the sheath until the tip of the dilator exits the sheath. The hub of the dilator may be locked into position.

Insertion Use and Removal of the Nagare™ Bi-Directional Steerable Sheath

- Follow standard access procedures to gain access into the vasculature using .032" x 180 cm J-tip guidewire.
- Flush the sheath side port and dilator lumens with sterile saline solution.
- Verify that the sheath is in a neutral (non-deflected) position. Advance the dilator and sheath together over the wire to the desired location under fluoroscopy.

Caution: If resistance is met when advancing or withdrawing any device, determine the cause and correct before continuing with the procedure.

- Remove the dilator when access to the desired vasculature is achieved.

Caution: To minimize potential for creating a vacuum in the sheath, remove the components and make catheter exchanges slowly. Always aspirate when removing the dilator or catheter.

- Aspirate until blood return is achieved and then flush with saline. It is recommended to flush through the side port.
- Attach a heparinized sterile saline flush/drip to the side port to reduce embolic risk or continue to flush through the side port while the sheath is in the vasculature.
- To insert a device through the Nagare™ Bi-Directional Steerable Sheath, pass the device through the hemostasis valve and advance it through the body of the sheath.

Caution: Do not manipulate the sheath in the heart without a device extending from the tip.

- To deflect the tip of the Nagare™ Bi-Directional Steerable Sheath, rotate the deflection knob. Note the position of the tip when utilizing the deflection mechanism.

Rotation of the deflection knob in a clockwise direction deflects the sheath tip towards the side port. Rotation of the deflection knob in a counter-clockwise direction deflects the sheath tip away from the side port. Perform all sheath manipulations under fluoroscopy.

- Upon completion of the procedure, ensure the steerable sheath is in the neutral position prior to removal. Follow standard hospital practices for vascular access closure.

Caution: Prior to withdrawal, always ensure the steerable sheath is in the neutral (non-deflected) position; this position is achieved by rotating the deflection knob in a direction to straighten the tip.

Glosary of Terms:

The symbols glossary table is located on the last page of this Instructions for Use.

DISCLAIMER OF WARRANTY

Note: Although the Nagare™ Bi-Directional Steerable Sheath, hereafter referred to as "product," has been manufactured under carefully controlled conditions, Kalila Medical and its affiliate, hereafter referred to as "Kalila," have no control over conditions under which this product is used. Kalila, therefore, disclaims all warranties, both expressed and implied with respect to the product, including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Kalila shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind Kalila to any representation or warranty with respect to the product. The exclusions and limitations set out above are not intended to and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected.

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