Glidesheath Slender® Tibial Pedal Kit

For lower extremity peripheral vasculature below the knee access
Hydrophilic Coated Introducer Sheath
(With Metallic Entry Needle)

INTRODUCER (With Metallic Entry Needle)
INTRODUCTEUR (Avec aiguille d'introduction métallique)
INTRODUCTOR (Aguja metálica de inserción)
INTRODUTOR (com agulha metálica de entrada)

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


MANUFACTURER: TERUMO MEDICAL CORPORATION
950 ELKTON BLVD., ELKTON MD, 21921, U.S.A

Registered Trademark M Coat and Glidesheath Slender are trademarks of TERUMO CORPORATION.
Marque Enregistrée M Coat et Glidesheath Slender sont des marques de TERUMO CORPORATION.
Marca Registrada M Coat y Glidesheath Slender son marcas comerciales registradas de TERUMO CORPORATION.
Marca Registrada M Coat e Glidesheath Slender são marcas comerciais da TERUMO CORPORATION.

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ENGLISH

Read all instructions before use.

INDICATION FOR USE
The Glidesheath Slender Tibial Pedal Kit is indicated to facilitate placing a catheter through the skin into the lower extremity peripheral vasculature below the knee.

CONTRAINDICATIONS
There are no known contraindications for the Glidesheath beyond standard diagnostic and/or interventional procedure contraindications.

POTENTIAL COMPLICATIONS
• Bleeding complications
• Distal embolization
• Infection and pain at puncture site
• Hematoma
• Intra-vessel thrombosis
• Arterial spasm
• Arterial injury/perforation/dissection
• Arterial embolism/occlusion
• Hemorrhage and hemorrhagic shock

DESCRIPTION/DESIGNATION
The Glidesheath Slender Tibial Pedal Kit consists of an introducer (a sheath and a dilator), a flex tip mini guide wire, and an entry needle. The sheath surface of this product is coated with "M Coat™", Terumo Corporation's exclusive hydrophilic coating. When the surface of the sheath is wet, it becomes very lubricious and the coefficient of friction is largely reduced.

< COMPONENT DESCRIPTION >
Refer to product label for appropriate system components.
The Glidesheath Slender is not made with natural rubber latex.
All components that enter blood vessels are radiopaque.
Sheath: Incorporates a 1-way valve and a 3-way stopcock connected by a side tube. The sheath can be used with a catheter of the same Fr. size or up to two Fr. sizes smaller without blood leakage at the 1-way valve.
Dilator: The precise fit of the dilator in the sheath allows for simultaneous motion of both dilator and sheath.

Mini Guide Wire
Entry Needle: Metallic Entry Needle
Guide Insertor

WARNING/PRECAUTIONS
WARNINGS
• Do not use a plastic guide wire with a metallic entry needle. Withdrawing the plastic wire through the metallic entry needle or advancing the entry needle over the plastic wire may cause the plastic part to shear, that may necessitate retrieval.
• Use of alcohol, antiseptic solutions or other solvents must be avoided, as they may adversely affect the surface of the sheath.
• Apply appropriate anticoagulant therapy to the patient.
• Do not reshape the tip of the mini guide wire.
• Manipulate the mini guide wire slowly and carefully not to damage the vessel wall, while monitoring the tip position and movement under image guidance.
PRECAUTIONS
• Before use, confirm the instructions for use of the drugs and devices to be used along with this device. This kit must be used by a trained physician under image guidance.
• Do not manipulate intraluminal devices if any resistance is met. Failure to exercise proper caution may lead to deformation of the sheath tube due to its thin wall, resulting in damage to the vessel.
• This product has been sterilized by ethylene oxide gas. For single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
• Remove the sheath from packaging carefully to avoid damage to the sheath tubing.
• Use care when handling or adjusting the sheath to avoid damage to the sheath tubing.
• Do not use if the unit package or the product has been damaged or soiled.
• Use immediately after opening the unit packaging.
• Dispose of safely after single use to avoid risk of infection.
• Before use, make sure the sheath (Fr.) and dilator size are appropriate for the access vessel and the system to be used.
• The entire procedure from skin incision to sheath removal must be carried out aseptically.
• The mini guide wire in this kit is designed to be used only for sheath insertion. Usage for any other purpose may result in damage to and/or separation of the wire, which may have to be removed from the vessel.

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• Do not use power injector for contrast media injection from the side tube while a device such as a dilator or a catheter is inside of the introducer sheath.
• Do not flush liquid such as contrast media or heparinized saline from side tube while a device such as a dilator or a catheter is inside of the introducer sheath.
• Do not inject drugs including oil components such as lipid emulsion, castor oil, interfacial active agent or solubilization agent such as alcohol, through the side tube. It may cause cracks on the stopcock.
• Be careful not to cut the side tube when holding it with forceps, or not to cut with scissors and knives.
• Do not excessively pull, push or bend the side tube and the stopcock.
• Do not pinch a side tube with forceps. It will cause scratches on it. Attention should be paid not to damage the cannula with forceps or sharp edged tools.
• Do not scratch the sheath with needle point, cutting tool, or other edged tools.
• Do not incline a guide wire and/or a catheter while inserting them through the introducer sheath valve.

INSTRUCTIONS FOR USE
1. If deemed necessary by the physician make a small skin incision at the puncture site with a surgical knife.
2. Insert the metallic entry needle into the artery (Fig. 1).
3. Before inserting mini guide wire, set the guide inserter at the metallic entry needle hub for easy insertion. Insert the selected flexible end of the mini guide wire through the metallic entry needle into the artery (Fig. 2).
CAUTION Do not withdraw the guide wire through the needle, as shearing of the guide wire may result. If resistance is met, do not advance or withdraw the mini guide wire until the cause of resistance is determined.
4. Remove the metallic entry needle over the mini guide wire.
CAUTION Dispose the metallic entry needle safely after use, being careful to avoid infections.
5. Connect a flushing line to the 3-way stopcock of the sheath. Fill the sheath assembly completely with heparinized saline, removing all air, also wet the surface of the sheath with heparinized saline solution.
CAUTION Do not turn the three way stop cock more than 180 degrees. The cock could misalign or come off and could cause blood leakage or it could shut down the path for the drug.
6. Prime the dilator using a syringe with heparinized saline.
7. Insert the vessel dilator fully into the sheath. The female hub of the sheath connects with the male hub of the dilator, and locks in place by means of grip (Fig. 3).
CAUTIONS
• Insert the dilator into the valve center of sheath. Forced dilator insertion missing the valve center may cause valve damage, resulting in blood leakage.
• Lock the dilator hub into the sheath hub securely. If the dilator hub is not locked into the sheath hub, only the sheath will advance into the vessel and the tip of the sheath may damage the vessel.
• Use sheath in a wet state. If the surface has dried, wet before use.
8. Insert the dilator and sheath together over the mini guide wire, and into the artery (Fig. 4).
9. "Unlock" the dilator hub from the sheath hub by bending the dilator hub downward (Fig. 5).
10. Slowly remove the dilator and mini guide wire, leaving the outer sheath in the artery (Fig. 6).
CAUTION If injection or sampling is necessary at this point, remove the mini guide wire, and use the dilator hub as an injection port before removing it. Slowly remove the dilator from the outer sheath. Rapid withdrawal of the dilator may result in incomplete closing of the 1-way valve, resulting in blood flow through the valve. If this occurs, replace dilator into the sheath and remove again slowly. After removing the dilator and mini guide wire, be careful while advancing the sheath. It may cause damage to the vessel.
11. Insert the catheter into the valve center of sheath.
12. Insert a catheter through the sheath and into the artery, and advance to the desired location (Fig. 7). When exchanging catheters, remove the used catheter and repeat step 12.
CAUTIONS
• Before removing or inserting the catheter through the sheath, aspirate blood from the 3-way stopcock to remove any fibrin deposition which may have accumulated in or on the tip of the sheath (Fig. 8).
• Be careful for rapid aspiration by a syringe. It may cause air to be drawn through the valve.
• When puncturing, suturing, or incising the tissue near the sheath, be careful not to damage the sheath.
• Do not clamp the sheath or bind with thread.
The flushing line may also be used as a continuous infusion site by connecting an infusion line to the 3-way stopcock (Fig. 9).
13. When inserting, manipulating or withdrawing a catheter from the sheath always hold the sheath in place.
CAUTION Do not place suture on the sheath tubing since this may restrict access/flow through the sheath.
14. After the intended procedure is completed, remove the catheter and then the sheath.

CAUTION FOR STORAGE
DO NOT STORE AT EXTREME TEMPERATURES AND HUMIDITY. AVOID DIRECT SUNLIGHT.

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SLIN0002 © TERUMO CORPORATION 2018-03-15