
Glidesheath™




Hydrophilic Coated Introducer Sheath
Gaine D'Introducer Enduit Hydrophile
Vaina Introducutora Revestido Hidrofílico
Bainha Introdutora Revestido Hydrophilic

MANUFACTURER: TERUMO MEDICAL CORPORATION
950 ELKTON BLVD, ELKTON MD, 21921, U.S.A

®: Registered Trademark M Coat and Glidesheath are trademarks of TERUMO CORPORATION.
®: Marque Enregistrée M Coat et Glidesheath sont des marques déposées d'ETERUMO CORPORATION.
®: Marca Registrada M Coat y Glidesheath son marcas registradas de TERUMO CORPORATION.
®: Marca Registrada M Coat e Glidesheath sao marcas registradas de TERUMO CORPORATION.

TERUMO

RFIN0075 © TERUMO CORPORATION 10 / 2015

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.
ATTENTION : la loi fédérale en vigueur aux États-Unis restreint la vente de ce dispositif uniquement par un médecin ou sur prescription médicale.
PRECAUCIÓN: Las leyes federales (EEUU) restringen la venta de este dispositivo por un médico o por orden de este.
ATENÇÃO: A lei federal (E.U.A.) restringe a venda deste dispositivo a médicos ou sob ordem médica.

PM-05028

ENGLISH

Read all instructions before use.

INDICATION FOR USE

The Glidesheath is indicated to facilitate placing a catheter through the skin into a vein or artery including but not limited to the radial artery.

CONTRAINDICATIONS:

There are no known contraindications for the Glidesheath beyond standard diagnostic and/or interventional procedure contraindications.

DESCRIPTION/DESIGNATION

The Glidesheath™ consists of an introducer (a sheath and a dilator), and a mini guide wire. Glidesheath™ products have an exclusive hydrophilic coating on the surface of the sheath. When the surface of the sheath is wet, it becomes very lubricous and the coefficient of friction is largely reduced. The Glidesheath is used to facilitate placing a catheter through the skin into a vein or artery including but not limited to the radial artery.

The Entry Needle is an accessory device which is used to gain access to the vein or artery, including but not limited to the radial artery, for placement of the mini guide wire.

The Mini Guide Wire is an accessory device which is used for placement of the sheath into the vein or artery including but not limited to the radial artery.

The Guide wire inserter is an accessory device which is used for placement of the mini guide wire into the entry needle.

< COMPONENT DESCRIPTION >

Refer to product labelling for appropriate system components.

All components that enter blood vessels are radiopaque

Sheath : Incorporates a 1-way valve and a 3-way stopcock connected by a side tube. This highly flexible sheath is designed to resist kinking.

Dilator : The precise fit of the dilator in the sheath allows for simultaneous motion of both dilator and sheath.

Mini Guide Wire : Material options include Stainless Steel, Plastic Jacketed, and Nitinol

Syringe : 2.5 ml flushing syringe Entry Needle : Metallic or Inner metallic needle and a plastic cannula

Guide Wire Inserter

WARNING/PRECAUTIONS

WARNINGS

- Do not use a plastic jacketed mini guide wire with a metallic entry needle. Withdrawing the plastic wire through the metallic entry needle or advancing the entry needle over the plastic jacketed wire may cause the plastic part to shear, that may necessitate retrieval.
- Use of alcohol, antiseptic solutions or other solvents must be avoided, as they may adversely affect the surface of the sheath.
- Apply appropriate anticoagulant therapy to the patient.
- Do not reshape the tip of the mini guide wire.
- Manipulate the mini guide wire slowly and carefully not to damage the vessel wall, while monitoring the tip position and movement under fluoroscopy.

PRECAUTIONS

- Before use, confirm the instructions for use of the drugs and devices to be used along with this device.
- This kit must be used by a trained physician under fluoroscopy.
- Prior to beginning radial artery access, the Allen test should be performed to assess the presence/adequacy of dual arterial circulation to the hand. Radial artery access is not recommended for patients with abnormal Allen test results.
- This kit is for single use only. Do not resterilize. Do not reuse. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- Do not use if the unit package or the product has been damaged or soiled.
- Use immediately after opening the unit packaging.
- Dispose of safely after single use to avoid risk of infection.
- Before use, make sure the sheath(Fr) and dilator size are appropriate for the access vessel and the system to be used.
- The entire procedure from skin incision to sheath removal must be carried out aseptically.
- Do not use power injector for contrast media injection from the side tube.
- Do not inject drugs including oil components such as lipid emulsion, castor oil, interfacial active agent or solubilization agent such as alcohol, through the side tube. It may cause cracks on the stopcock.
- When the catheter is inserted in the sheath, do not flush from the 3-way stopcock.
- Be careful not to cut the side tube when holding it with forceps, or not to cut with scissors and knives
- Do not excessively pull, push or bend the side tube and the stopcock.

CAUTION : Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

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INSTRUCTIONS FOR USE

- If deemed necessary by the physician make a small skin incision at the puncture site with a surgical knife.
- Insert an entry needle to puncture the blood vessel (Fig. 1A and 1B). When using a metallic entry needle with plastic cannula remove inner metallic needle before wire insertion (Fig. 1A).
CAUTION Do not reinsert the inner metallic needle into the plastic cannula. Dispose of the removed inner metallic needle safely being careful to avoid infections.
- Insert the flexible end of the mini guide wire through the plastic cannula into the vessel (Fig. 2A and 2B). If using a guide wire inserter set it at the entry device hub (metallic needle or plastic cannula) for easy insertion.
CAUTION Advance or withdraw the mini guide wire slowly. If resistance is met, do not advance or withdraw the mini guide wire until the cause of resistance is determined.
- Remove the entry device over the mini guide wire.
CAUTION Dispose of the removed entry device safely being careful to avoid infections.
- Connect a flushing line to the 3-way stopcock of the introducer sheath. Fill the sheath assembly completely with heparinized saline, removing all air. Wet the surface of the sheath with heparinized saline solution.
CAUTION Do not turn the three way stopcock more than 180 degrees. The stopcock could misalign or come off and could cause blood leakage or it could shut down the path for the drug.
- Prime the dilator using the syringe with heparinized saline.
- Insert the vessel dilator fully into the sheath. The female hub of the sheath connects with the male hub of the dilator, and laces in place by means of grip (Fig. 3).
CAUTIONS
 - Insert the dilator into the valve center of sheath. Forced dilator insertion missing the valve center may cause valve damage, resulting in blood leakage.
 - Lock the dilator hub into the sheath hub securely. If the dilator hub is not locked into the sheath hub, only the sheath will advance into the vessel and the tip of the sheath may damage the vessel.
 - Use Glidesheath™ introducer sheaths in a wet state. If the surface has dried, wet before use.
- Insert the dilator and sheath together over the mini guide wire, and into the blood vessel (Fig. 4).
- "Unlock" the dilator hub from the sheath hub by bending the dilator hub downward (Fig. 5).
- Slowly remove the dilator and mini guide wire, leaving the outer sheath in the vessel (Fig. 6).
CAUTION If injection or sampling is necessary at this point, remove the mini guide wire, and use the dilator hub as an injection port before removing it. Slowly remove the dilator from the outer sheath. Rapid withdrawal of the dilator may result in incomplete closing of the 1-way valve, resulting in blood flow through the valve. If this occurs, replace the dilator into the sheath and remove again slowly. After removing the dilator and mini guide wire, be careful for advancing the sheath. It may cause damage to the vessel.
- Insert the catheter into the valve center of sheath.
CAUTIONS
 - Before inserting the catheter through the sheath, aspirate blood from the 3-way stopcock to remove any fibrin deposition which may have accumulated in or on the tip of the sheath.
- Insert a catheter through the sheath and into the blood vessel, and advance to the desired location (Fig. 7). When exchanging catheters, remove the used catheter and repeat step12.
CAUTIONS
 - Before removing or inserting the catheter through the sheath, aspirate blood from the 3-way stopcock to remove any fibrin deposition which may have accumulated in or on the tip of the sheath (Fig. 8).
 - Be careful for rapid aspiration by a syringe. It may cause air to be drawn through the valve.
 - When puncturing, suturing, or incising the tissue near the sheath, be careful not to damage the sheath. Do not put a clamp on the sheath nor bind it with a thread.
- The flushing line may also be used as a continuous infusion site by connecting an infusion line to the 3-way stopcock (Fig. 9).
- When inserting, manipulating or withdrawing a catheter from the sheath always hold it in place. To temporarily suture the sheath (for continued access) use the suture eye.
CAUTION Do not place suture on the sheath tubing since this may restrict access / flow through the sheath.
- After the intended procedure is completed, remove the catheter and then the sheath.

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