GLIDEWIRE® GT

ENGLISH

Please read all instructions prior to use.

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

INDICATIONS

The GLIDEWIRE GT is designed to direct a catheter to the desired anatomical location in the peripheral vasculature during diagnostic or interventional procedures. This device is not intended for neurovascular or coronary interventions.

DESCRIPTION

The GLIDEWIRE GT consists of a core wire of a superelastic alloy coated with a hydrophilic polymer material and has a flexible radiopaque gold coil around the distal tip of the core wire. The RADIFOCUS® TORQUE DEVICE supplied facilitates guide wire manipulation.

CONTRAINDICATIONS

Do not use the GLIDEWIRE GT in percutaneous transluminal coronary angioplasty (PTCA) or percutaneous transluminal angioplasty (PTA) for severe stenosis. Advancing or withdrawing the wire against large resistance may result in torsion or bending of the wire's tip and subsequently in damage to the vessel and/or damage to or separation of the wire.

COMPLICATIONS

Possible complications associated with the use of the GLIDEWIRE GT include, but are not limited to, the following :

·vessel perforation ·intimal tear ·vessel damage ·vessel spasm ·hematoma at the puncture site ·infections

WARNINGS/PRECAUTIONS

Failure to abide by the following warnings might result in damage to the vessel, abrasion of the hydrophilic coating, release of plastic fragments from the GLIDEWIRE GT, damage to or breakage/separation of the wire, that may necessitate retrieval.

- WARNINGS
- Do not manipulate and/or withdraw the GLIDEWIRE GT through a metal entry needle, a metal dilator, or a metal guide wire inserter. Manipulation and/or withdrawal through a metal entry needle, a metal dilator, or a metal guide wire inserter may result in destruction and /or separation of the outer polyurethane coating requiring retrieval. A plastic entry needle is recommended when using this wire for initial placement.
 - Do not use the GLIDEWIRE GT with a metal-tipped catheter, a laser catheter, or a catheter that has a rotating cutting design such as an atherectomy catheter. Use of the GLIDEWIRE GT with such catheters may result in damage to the surface of the wire and/or separation of the wire.
 - Do not attempt to reshape the GLIDEWIRE GT unless it is indicated as "shapeable" on the package. Attempting to reshape the wire may cause damage, resulting in the release of wire fragments into the vessel.
 - When reshaping the shapeable GLIDEWIRE GT tip use the mandrel. Do not use heat, sharp metal instruments, or excessive force. These methods may result in abrasion of hydrophilic coating and damage to the wire.
 - A retrieving device, such as a gripper or basket forceps, can only be used after the GLIDEWIRE GT has been removed from the patient's vessel. Using a retrieving device while the GLIDEWIRE GT is in the vessel may cause the GLIDEWIRE GT to break.
 - Manipulate the GLIDEWIRE GT slowly and carefully in the vessel while confirming the behavior and location of the wire's tip under fluoroscopy to avoid damage to the vessel. If any resistance is felt, the behavior and/or position of the wire's tip seems improper, the wire is kinked, or vessel spasm is suspected, stop manipulating the wire (and the catheter) and determine the cause carefully by fluoroscopy to take suitable remedial actions. Then remove the wire slowly without turning. Failure to exercise proper caution may result in bending, kinking, separation of the guide wire's tip, damage to the catheter, or damage to the vessel.
 - Do not turn the proximal end of the GLIDEWIRE GT three or more turns in succession in the same direction if the wire's tip is trapped due to vessel spasm or some other cause. Separation of the wire may result.
 - Carefully manipulate the guiding catheter and observe it for undesirable behaviors while the GLIDEWIRE GT is in the vessel. Especially when the wire's tip is trapped due to vessel spasm or some other cause, damage to or separation of the wire and/ or damage to the vessel may result if the guiding catheter slips out of position and exerts too large a load on the wire.
 - Do not attempt to use the GLIDEWIRE GT if it has been bent, kinked or damaged. Use of a damaged wire may result in damage to the vessel or the release of wire fragments into the vessel.
 - The GLIDEWIRE GT may slide entirely into the catheter or slip out of the catheter because of its low sliding friction.
 - The GLIDEWIRE GT should be used by a physician, who is well trained in manipulation and observation of guide wires under fluoroscopy.
 - The GLIDEWIRE GT has been sterilized by ethylene oxide gas. For single use only. **Do not** resterilize or reuse.
 - Sterile, non-toxic and non-pyrogenic in unopened and undamaged unit package. **Do not** use if the unit package or the guide wire is broken or soiled.
 - The GLIDEWIRE GT should be used immediately after opening the package and be disposed of safely and properly after use, following local regulations for medical waste management.
 - The entire operation should be carried out aseptically.
 - When using a drug or a device concurrently with the GLIDEWIRE GT, the operator should have a full understanding of the properties/ characteristics of the drug or device so as to avoid damage to the GLIDEWIRE GT.
 - The surface of the GLIDEWIRE GT **is not** lubricious unless it is wet. Before taking it out of its holder and inserting it through a catheter, fill the holder and the catheter with heparinized physiological saline solution.
 - When reinserting the GLIDEWIRE GT back into the holder, take care **not** to damage the wire's hydrophilic polymer coating with the edge of the holder.
 - Use the enclosed RADIFOCUS TORQUE DEVICE to handle the GLIDEWIRE GT. **Do not** use a metal torque device with the GLIDEWIRE GT. Use of a metal torque device may result in damage to the wire. Also, **do not** slip a tightened up torque device over the wire, as this may result in damage to the wire.
 - Due to variations of certain catheter tip inner diameters, abrasion of the hydrophilic coating may occur during manipulation. If any resistance is felt during introduction of the catheter, it is advisable to stop using such catheters.
 - **Do not** manipulate the GLIDEWIRE GT through a tightened up rotating hemostasis valve, as this may result in damage to the wire.
 - After removal from the patient's vessel, and prior to reinserting it into the same patient during the same catheterization, the GLIDEWIRE GT should be rinsed in a bowl full of heparinized physiological saline solution. Any blood residues still adhering to the wire can be removed by wiping once with a gauze moistened with heparinized physiological saline solution. Use of alcohol, antiseptic solutions or other solvents must be avoided, because they may adversely affect the surface of the GLIDEWIRE GT.

DIRECTIONS FOR USE

1. Remove the GLIDEWIRE GT and the holder together from the tray.
2. Fill the holder with heparinized physiological saline solution through the hub of the holder using a syringe.
3. Remove the inserter and then the GLIDEWIRE GT from the holder, and inspect the wire prior to use to verify it is lubricated. If the GLIDEWIRE GT **can not** be easily removed from the holder, inject more heparinized physiological saline solution into the holder and try again.
4. Prime the catheter with heparinized physiological saline solution to ensure smooth movement of the GLIDEWIRE GT within the catheter.
5. Keep at least 5 cm of the wire extended out of the hub of the catheter during introduction. A RADIFOCUS TORQUE DEVICE, which is included in the package, is recommended for easier handling/manipulation of the wire.
6. When using the Torque Device, loosen the nut and slip it over the proximal end of the GLIDEWIRE GT.
7. When the Torque Device is in the desired position on the GLIDEWIRE GT, tighten the nut. The Torque Device may be repositioned or removed by loosening the nut and retightening it.
8. Manipulate the GLIDEWIRE GT in the vessel while monitoring the behavior and location of the wire's tip under fluoroscopy.

- WARNING
- If the tip of the GLIDEWIRE GT should be trapped, **do not** turn or draw it by force. Instead determine the cause to take suitable remedial actions, then withdraw the wire slowly without turning. Failure to exercise proper caution may result in damage to the wire's tip.

PRECAUTION FOR STORAGE

Avoid exposure to water, direct sunlight, extreme temperatures and high humidity during storage. Store under controlled room temperature.

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TERUMO is a trademark of TERUMO CORPORATON.
® : Registered Trademark. ® : Marque Enregistrée. © TERUMO CORPORATION 2019 - 11