TR Band®

Radial Artery Compression Device
 Dispositif de compression de l'artère radiale
 Kompressionsvorrichtung für die Radialarterie
 Dispositivo de compresión arterial radial
 Dispositivo de Compressão da Artéria Radial
 Dispositivo di compression dell'arteria radiale
 Compressie-instrument voor radiale slagader
 Radiell artärkompressionsenhet.
 Udstyr til kompression af arteria radialis
 Komprimeringsenhet for radial arterie
 Värntinåvaltimon puristuslaite
 Συσκευή συμπίεσης για κερκιδική αρτηρία
 Устройство для гемостаза лучевой артерии
 Urządzenie do ucisku tętnicy promieniowej

Radiális artériás kompressziós eszköz
 Kompresni prostředek pro a. radialis
 Kompresná pomůcka na arteria radialis
 Radyal Arter Kompresyon Cihazı
 Kodarluuarteri kompressioonisead
 Radiālās artērijas kompresijas ierīce
 Radialinės arterijos suspaudimo įrenginys
 Pripomoček za kompresijo radialne arterije
 Uređaj za kompresiju radialne arterije
 Dispozitiv de compresie a arterei radiale
 Пристрій для компресії променевої артерії
 Uređaj za kompresiju radialne arterije
 Изделие за компресия на радиална артерия

®: Registered Trademark
 ®: marque déposée
 ®: Eingetragenes Markenzeichen
 ®: Marca comercial registrada
 ®: Marca comercial registada
 ®: Marchio registrato
 ®: Geregistreerd handelsmerk
 ®: Registrerat varumärke
 ®: registreret varemærke
 ®: Registrert varemerke
 ®: Rekisteröity tavaramerkki
 ®: Σημειωμένο
 ®: Зарегистрированный товарный знак
 ®: Zastrzeżony znak towarowy
 ®: Bejegyzett védjegy
 ®: Registrovaná ochranná známka
 ®: Registrovaná ochranná známka
 ®: Tescilli Ticari Marka
 ®: Registreeritud kaubamärk
 ®: Registrētā preču zīme
 ®: Registruotais preču zīmējums
 ®: registrirana blagovna znamka
 ®: Registrovana robna marka
 ®: Marcă înregistrată
 ®: Регистрирана търговска марка
 ®: Зареєстрована торгова марка
 ®: Registrirani zaštitni znak

TR Band is a trademark of TERUMO CORPORATION.
 TR Band est une marque commerciale de TERUMO CORPORATION.
 TR Band ist ein eingetragenes Markenzeichen von TERUMO CORPORATION.
 TR Band es una marca registrada de TERUMO CORPORATION.
 TR Band é uma marca comercial da TERUMO CORPORATION.
 TR BAND è un marchio commerciale di TERUMO CORPORATION.
 TR BAND is een handelsmerk van TERUMO CORPORATION.
 TR band är ett varumärke som tillhör TERUMO CORPORATION.
 TR Band er et varemærke tilhørende TERUMO CORPORATION.
 TR Band er et varemerke for TERUMO CORPORATION.
 TR Band on TERUMO CORPORATION -yhtiön tavaramerkki.
 Το TR Band είναι εμπορικό σήμα της TERUMO CORPORATION.
 TR Band является торговой маркой TERUMO CORPORATION.
 TR Band jest znakiem towarowym firmy TERUMO
 A TR Band a TERUMO CORPORATION védjegye.
 TR Band is ochranná známka společnosti TERUMO CORPORATION.
 Značka TR Brand je ochrannou známkou společnosti TERUMO CORPORATION.
 TR Band, TERUMO CORPORATION in ticari markasıdır.
 TR Band on ettevõtte TERUMO CORPORATION kaubamärk.
 TR Band ir TERUMO CORPORATION preču zīme.
 „TR Band” yra TERUMO CORPORATION prekės ženklas.
 Pas TR Band je blagovna znamka družbe TERUMO CORPORATION.
 TR Band je robna marka kompanije TERUMO CORPORATION.
 TR Band este o marcă comercială a TERUMO CORPORATION.
 TR Band е търговска марка на TERUMO CORPORATION.
 TR Band е търговска марка TERUMO CORPORATION.
 TR Band は TERUMO CORPORATION の商標です。
 TR band är ett varumärke som tillhör TERUMO CORPORATION

TERUMO MEDICAL CORPORATION

950 Elkton Blvd., Elkton,
 Maryland 21921, U.S.A

MADE IN U.S.A.

Green marker: Align this marker to the puncture site
 Port d'injection d'air avec vanne
 Lufteinblasöffnung mit Ventil
 Puerto de inyección de aire con válvula
 omada de Injeção de Ar com Válvula
 Porta di iniezione dell'aria con valvola
 Luchtinjectieopening met klep
 Luftinjektionsport med ventil
 Lufttilførselsport med ventil
 Luftinjektionsport med ventil
 Venttiiläilä varustettu ilmanruiskutusportti
 Όρος έγχυσης αέρα με βαλβίδα
 Порт с клапаном для подачи воздуха
 Port pompowania powietrza z zaworem
 Levegőbefecskendező port szeleppel
 Port na vstříkávání vzduchu s ventilem
 Port na aplikáciu vzduchu s ventilom
 Valfli Hava Enjeksiyon Portu
 Klapięga õhusustimisport
 Gaisa iesmidzināšanas pieslēgvieta ar vārstu
 Oro švirkštimo anga su vožuvu
 Odrpina za injiciranje zraka z ventnilom
 Otvor za ubrizgavanje vazduha sa ventnilom
 Port de injectare a aerului cu supară
 Отвір для впрскування повітря з клапаном
 Otvor za ubrizgavanje zraka s ventnilom
 Otvor za injektiranje na vъздyх с клапан

Inflator Tube Clip
 Tube de gonflage
 Schlauch für Aufblasvorrichtung
 Tubo del inflador
 Tubo de Insuflação
 Tubo del dispositivo di gonfiaggio
 Pomsilang
 Uppblåsningsslang
 Klemme til oppustningsrør
 Oppumpingsrør
 Täytöputki
 Σωληνιοκόξ διόγκωσής
 Трубка насоса
 Rurka pompki
 Felfújoeszköz csőve
 Hadička inflátoru
 Hadička inflátora
 Sijmre Hortumu
 Taitmisvoolik
 Pēpõšanas ierices caurulite
 Pipitõtkilo vamzdelis
 Cevka za inficacio
 Tub de umflare
 Трубка для подавання повітря
 Сјев параве за напуњавање
 Тръбичка за раздуване


Air injection port with valve
 Marqueur vert (aligner ce marqueur sur le site de ponction)
 Grüner Marker (richten Sie diesen Marker auf die Einstichstelle aus)
 Marcador verde (alineare este marcador com o local da punção)
 Marcador Verde (alinhar este marcador com o local da punção)
 Marker verde (allineare questo marker al sito di puntura)
 Groene markering (breng deze markering op één lijn met de prikplek)
 Grön markör (passa in denna markör med punktionsstället)
 Grøn marker (placér denne marker overfor punkturstedet)
 Grønn marker (juster denne markøren til punkteringsstedet)
 Vihreä merkki (kohdistä tämä merkki pistokohtaan)
 Πρόσόνος δείκτης (επιθυραμίστε αυτόν τον δείκτη με το σημείο της ποράκέντησης)
 Зелёная метка (соотнесите эту метку с местом прокола)
 Zielony znacznik (wyrównać ten znacznik z miejscem wkłucia)
 Zöld jelölés (igazítsa a jelölést a szúrás helyéhez)
 Zelená značka (zarovnejte tuto značku k místu punkcie)
 Zelená značka (zarovnajte túto značku s miestom punkcie)
 Yeşil işaret (bu işareti ponksiyon bölgesiyle hizalayın)
 Roheline marker (joondage see marker punktisoonikohaga)
 Žaļa atzīme (izlīdziniet šo atzīmi ar punkcijas vietu)
 Žalias žymeklis (sulygiuokite šį žymeklį su pradūrimo vieta)
 Zelena oznaka (to oznako poravnajte z mestom vboda)
 Zeleni marker (poravnajte ovaj marker sa mestom uboda)
 Marčaj verde (aliniati acest marčaj la locul punctiei)
 Зелений маркер (що наводиться на місце проколу)
 Zelena oznaka (poravnajte ovu oznaku s mjestom uboda)
 Зелен маркер (подравняйте този маркер с мястото на пункция)

Adjustable fastener
 Fixation réglable
 Einstellbarer Verschluss
 Sujeción ajustable
 Fecho Ajustável
 Chiusura regolabile
 Versteerbare sluiting
 Luftinjektionsport med ventil
 Indstillelig fastgørelse
 Justerbar feste
 Säädettävä kiinnike
 Ρυθμιζόμενο κούμπωμα
 Регулируемая застёжка
 Regulowane zapięcie
 Állítható rögzítő
 Nastavitelný fixátor
 Nastavitelný fixátor
 Ayarlanabilir Tutturuç
 Reguleeritav kinniti
 Regulējams stiprinājums
 Reguliojamas užsegiklis
 Nastavljiva zaponka
 Podesiivi pričvršćivač
 Element de fixare reglabil
 Регульована застібка
 Podesiivi zatvarač
 Регулируема закопчалка

Compression balloon (small)
 Ballon de compression (petit)
 Kompressionsballon (klein)
 2Globo de compresión (pequeño)
 Balão de Compressão (pequeno)
 Palloncino a compressione (piccolo)
 Compressieballon (klein)
 Kompressionsballon (liten)
 Kompressionsballon (lille)
 Komprimeringsboble (liten)
 Puristusilmapallo (pieni)
 Μπαλόνι συμπίεσης (μικρό)
 Компрессионный баллон (маленький)
 Balon uciskowy (mały)
 Kompressziós ballon (kicsi)
 Kompresni balónek (malý)
 Kompresny balónik (malý)
 Kompresyon Balonu (küçük)
 2Kompressieballon (våike)
 Kompressijas balons (mazs)
 Suspaudimo balionas (mažas)
 Kompresijski balon (velik)
 Kompresioni balon (mal)
 Balon de compresie (mic)
 Компресійна камера (маленька)
 Kompresijski balon (mal)
 Компресиращ балон (малък)

Compression balloon (large)
 Ballon de compression (grand)
 Kompressionsballon (groß)
 Globo de compresión (grande)
 Balão de Compressão (grande)
 Palloncino a compressione (grande)
 Compressieballon (groot)
 Kompressionsballon (stor)
 Kompressionsballon (stor)
 Komprimeringsboble (stor)
 Puristusilmapallo (iso)
 Μπαλόνι συμπίεσης (μεγάλο)
 Компрессионный баллон (большой)
 Balon uciskowy (duży)
 Kompressziós ballon (nagy)
 Kompresni balónek (veľký)
 Kompresny balónik (veľký)
 Kompresyon Balonu (büyük)
 Kompressieballon (suur)
 Kompressijas balons (liels)
 Suspaudimo balionas (didelis)
 Kompresijski balon (majhen)
 Kompresioni balon (velik)
 Balon de compresie (mare)
 Компресійна камера (велика)
 Kompresijski balon (velik)
 Компресиращ балон (голям)

Support plate
 Plaque de support
 Stützplatte
 Placa de soporte
 Placa de Apoio
 Piastra di supporto
 Steunplaat
 Stödplatta
 Støtteplate
 Støtteplate
 Tukilevy
 Πλάκα υποστήριξης
 Пластина-держатель
 Plytka podparowa
 Tartóleméz
 Podpurná diaha
 Podporná plocha
 Destek Plakası
 Tugiilaat
 Atbalsta plāksne
 Atraminė plokštė
 Podporna ploščica
 Potporna ploča
 Placă suport
 Опорна пластина
 Potporna pločica
 Опорна пластина

Belt	Adjustable fastener
Ceinture	Fixation réglable
Gürtel	Einstellbarer Verschluss
Correa	Sujeción ajustable
Fita	Fecho Ajustável
Cinghia	Chiusura regolabile
Riem	Verstelbare sluiting
Rem	Justerbar fästaneländning
Rem	Indstillelig fastgørelse
Belt	Justerbar feste
Hihna	Säädettävä kiinnike
Ιμνόςας	Ρυθμιζόμενο κούμπωμα
Ремешок	Регулируемая застёжка
Pasek	Regulowane zapięcie
Szj	Állítható rögzítő
Pás	Nastavitelný fixátor
Pás	Nastavitelný fixátor
Kayış	Ayarlanabilir Tutturucu
Riim	Reguleeritav kinniti
Josta	Regulējams stiprinājums
Dirželis	Reguliuojamas užsegiklis
Trak	Nastavljiva zaponka
Pojas	Podesivi pričvršćivač
Centură	Element de fixare reglabil
Ремінець	Регульована застібка
Pojas	Podesivi zatvarač
Колан	Регулируема закопчалка

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Radial Artery Compression Device

INDICATION FOR USE

The TR Band® is a compression device indicated to apply compression in order to achieve hemostasis while allowing the user to maintain patency of the radial artery after a Transradial procedure (patient hemostasis).

WARNINGS

- Do not inject air into any port other than the air injection port of this device. This could cause arterial or venous embolization.
- Do not inject more than 18mL of air. There is a possibility of damaging the device if this volume is exceeded.
- Do not soak or wipe with agents containing organic solvents. Such agents may cause damage to the TR Band.
- Patient should not be left unattended while the TR Band is in use.
- Do not use TR Band in MRI field.

PRECAUTIONS

(Rx ONLY) CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

- Do not use the TR Band Inflator for other purposes than to inflate the TR Band. It is not designed for other purposes.
- The TR Band can be inflated only by using the TR Band Inflator included in the kit. If other devices are used, adequate air compression cannot be achieved.
- Only use air for inflation. Do not use other media for inflation.
- Depending on the patient's condition and the degree of balloon pressure, an adverse event including artery occlusion, hypodermic hematoma, bleeding, pain, or numbness may occur. Check the progress of the hemostasis and adjust the air pressure accordingly.
- When connecting the TR Band Inflator included in the kit to the TR Band, keep the plunger in place. Releasing the plunger will cause the air to expel from the TR Band. Loss of air compression may cause bleeding.
- Do not put excessive pressure on the tip of the TR Band Inflator; this may result in damage of the tip.
- If the inflator tip is broken while connected with the TR Band, it could cause air leakage and bleeding.
- This product has been sterilized by ethylene oxide gas. For single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- While using, be careful not to put external pressure on the device that could cause damage.
- Be careful that no foreign particles get into the air injection port when injecting air. The air could leak.
- If there is itching or redness of the skin while compressing, stop using and treat appropriately.
- At extremely low temperature, there is a possibility of damage due to decrease in resistance to shock.
- Take care not to damage the TR Band with needle, scissors and other sharp instrument. This may result in air leakage and bleeding.
- Do not use this device if the unit package or the product has been damaged or soiled.
- The device should be used immediately after opening the package and disposed of safely and properly after use.
- Care should be taken to prevent infection during manipulation.
- Avoid exposure to water, direct sunlight, extreme temperature, or high humidity during storage.

DIRECTIONS FOR USE

1. After the procedure, withdraw the sheath by 2-3 cm (Fig. 1).
2. Align the green marker, which is located on the center of the compression balloon (large), 1-2 mm proximal to the puncture site and fix the belt on the wrist with the adjustable fastener (Fig. 2).

CAUTIONS

- This device must be positioned differently when used on the left or right wrist. When attaching the device, ensure that the orientation of the palm symbol below the Terumo logo on the support plate matches patient's palm. (Fig. 2).
- Make sure the fastener is stable and not slanted.

Note: Use of pulse oximeter is suggested on the index finger or the thumb to record an adequate waveform.

3. To inflate the compression balloon, inject the appropriate volume of air using the TR Band Inflator, which is included in the kit (Fig. 3).

Nominal air injection volume: 13mL

Maximum air injection volume: 18mL.

WARNING

Prior to injection, confirm that you are injecting through the tube marked with this mark  and do not inject air through any other port.

4. Remove the sheath and confirm that there is no bleeding from the puncture site (Fig. 4). If bleeding is observed, inject more air (not exceeding a total of 18mL) until it stops.

Note: While manually occluding the ulnar artery, cessation of the oximetry waveform indicates complete compression of the radial artery. If the plethysmographic signal is absent, remove 1-2 mL from the compression bladder while maintaining hemostasis and manually occlude the ulnar artery until a pulsatile plethysmographic signal is observed. A pulsatile oximeter signal will be observed if radial artery patency has been achieved.

5. Check the progress of hemostasis and adjust the air pressure of the balloon with the TR Band Inflator accordingly (Fig. 5).

CAUTION

Bleeding could occur if the adjustable fastener comes off during compression. Depending on the patient's wrist size and the way the device is applied, the adjustable fastener could come off. Ensure the fastener is fixed properly to minimize risk of losing and loss of applied pressure which can cause bleeding.

6. Before removing, confirm that bleeding has stopped.

*Air injection volume and compression time may differ according to the patient's condition, the heparin volume, and the size of the puncture site. Check the puncture site and adjust accordingly.

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