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## DEVICE DESCRIPTION

The FineCross M3 consists of a catheter shaft, a hub and an anti-kink protector. The catheter shaft is comprised of an outer layer and an inner layer reinforced by braided mesh, and has a radiopaque marker at the distal tip. The distal part of the catheter outer surface is coated with a hydrophilic polymer, which becomes lubricious when wet.

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
# ENGLISH

Please read all instructions prior to use.

## INDICATIONS FOR USE

The product (FineCross M3) is intended to be percutaneously introduced into blood vessels and support a guide wire while performing PCI (percutaneous coronary intervention). The product is also intended for injection of radiopaque contrast media for angiography. The product should not be used in cerebral and peripheral vessels.

## CONTRAINDICATIONS

- Patients who had previous coronary artery spasm due to the possibility of acute coronary occlusion.
- Pregnancy or suspected pregnancy. X-ray exposure could damage fetus.
- Patients in whom CABG (coronary artery bypass graft) is not applicable. Emergent CABG is required in case of acute ischemic complications.
- This device should not be used in patients who cannot be systemically anticoagulated.

## COMPLICATIONS

Possible complications of PCI include, but are not limited to the following:

- Septicemia
- Cardiopulmonary arrest
- Fetus exposure
- Acute myocardial infarction
- Embolization
- Thrombosis
- Bleeding complications
- Perforation
- Coronary artery dissection
- Coronary artery spasm

## WARNINGS AND PRECAUTIONS

### WARNINGS

- **Carefully handle the product under high resolution fluoroscopy. If any resistance is felt while handling the product, immediately stop the manipulation and find out the cause to avoid damage to blood vessels and separation or breakage of the product.**
- **Before inserting / withdrawing the product, clean the surface of the guide wire with gauze moistened with saline solution. Advancing / withdrawing the product over a guide wire with residual blood on its surface or a guide wire which is not fully wet may result in separation or breakage of the product.**
- **Do not torque the product excessively while the distal part of the product crosses the stenosis or is in the stent.**
- **Be careful when inserting/withdrawing the product through an opening of the stent struts to avoid damage to the product. A scratch by the stent struts may result in separation or breakage of the product.**
- **Take extra care when exchanging the guide wires leaving the product in the arteries. Carefully insert a guide wire into the product. If any resistance is felt, stop the manipulation and remove the product together with the guide wire to avoid separation or breakage of the product.**
- **Perform appropriate anticoagulant or antiplatelet therapy according to the patient's condition in order to avoid complications, such as thrombotic embolization.**
- **Do not reshape the product by any means. Attempting to reshape the catheter may cause damage, resulting in the release of catheter fragments into the vessel.**
- **Failure to abide by the warnings in this labeling might result in damage to the device coating, which may necessitate intervention or result in serious adverse events.**
- **Do not use a power injector for the purpose of injection of contrast media.**

**PRECAUTIONS**

- This product has been sterilized by ethylene oxide gas. For single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- Do not apply agents containing organic solvents or oleaginous contrast media in order to avoid breakage of the product.
- Do not use the product for lesions in the left main trunk for which no compensation of blood flow by bypass or collateral circulation is available in order to avoid acute coronary occlusion.
- The product must be used by physicians who are well trained in PCI procedures.
- Do not use the product at institutions where an emergent CABG can not be performed in case of severe complications.
- Do not torque the product excessively if it is bent in order to avoid separation or breakage of the product.
- Do not soak the product in agents containing organic solvents, such as alcohol for disinfection. Do not clean the product with such agents. Failure to observe this precaution could damage or break the product or cause loss of lubricity.
- The entire procedure should be carried out aseptically.
- Carefully inspect the unit package to make sure that no damage to the sterile barrier has occurred.
- Sterile and non-pyrogenic in an unopened and undamaged unit package. Do not use if the unit package or the product has been damaged or soiled.
- The product should be used immediately after opening the package, and be disposed of safely and properly after use following the local regulations. The product is biohazardous as it is contaminated with blood.
- Maximum injection pressure should not exceed 2068 kPa (300psi).
- Do not use under MRI.
- Be careful when manipulating, advancing and/or withdrawing these devices through stents or other devices with sharp edges, or through tortuous or calcified blood vessels. These may result in destruction and/or separation of the outer coating, which may lead to clinical adverse events, resulting in coating material remaining in the vasculature or device damage.
- Avoid wiping the device with dry gauze as this may damage the device coating.
- Avoid excessive wiping of the coated device.
- Avoid using alcohol, antiseptic solutions, or other solvents to pre-treat the device because this may cause unpredictable changes in the coating which could affect the device safety and performance.
- Infusion flow rates for this product at maximum injection pressures for saline and iodixanol (11.2 mPa s 37°C) are provided below. Do not exceed the maximum permissible injection pressure. Contrast media should be injected at 37°C.

Usable Length	Flow Rate		Maximum Permissible Pressure
	Saline	Iodixanol	
130cm	1.48 mL / sec	0.18 mL / sec	300 psi (2,068 kPa)
150cm	1.41 mL / sec	0.16 mL / sec	300 psi (2,068 kPa)

## DIRECTIONS FOR USE

### 1. Preparation

- 1-1 Carefully remove the product in its holder from the package.
- 1-2 Soak the product in the holder in a heparinized saline solution bath.
- 1-3 Fill the holder with heparinized saline solution through the hub of the holder using a syringe, to thoroughly wet the surface of the product.

**CAUTION** The heparinized saline solution should be injected slowly into the holder so that the product is not driven out of its holder.

- 1-4 Carefully take the product out of the holder.

**CAUTION** Do not bend the product at the edge of the holder. The product may break or separate.

### 2. Insertion of the product

- 2-1 Prime the product with heparinized saline solution through the catheter hub in order to remove air inside.
- 2-2 Insert a guide wire (O.D. 0.014" (0.36 mm) or smaller) into the catheter hub and advance tip of the wire to the end of the catheter.

**CAUTION** Take care not to damage the product, when a guide wire is inserted from its distal end.

- 2-3 Open the hemostasis valve of the Y connector attached to the guiding catheter and insert the product into the valve.

**CAUTION** Make sure that the hemostasis valve is open enough for insertion of the product. If not, the valve may cause resistance.

- 2-4 Advance the distal end of the product to 2-3cm proximal to the end of the guiding catheter under high resolution fluoroscopy.

**CAUTION** Carefully insert the product into the guiding catheter with other devices in it. The other devices may be forced into the guiding catheter, which may result in damaging blood vessels.

- 2-5 Advance the guide wire into the target coronary artery and place the distal tip of the wire to the extreme periphery of the coronary artery under high resolution fluoroscopy. Inject a radiopaque dye in order to make sure that the guide wire crosses the stenosis. Perform the angiogram from multiple angles in order to make sure that the guide wire is placed at the target vessel.
- 2-6 Securely fix position of the guide wire and the guiding catheter and slowly advance the distal end of the product bit at a time over the guide wire until the radiopaque marker crosses the stenosis.

**CAUTIONS**

- Carefully handle the product inside coronary arteries because of its hydrophilic coating.
- Do not insert the product into small vessels, which have smaller inner diameter than outer diameter of the product.
- The product can be inserted into severely tortuous vessels but only at the physician's discretion.
- Do not torque the product if it is or seems stuck in order to avoid separation or breakage of the product.

- 2-7 In case injection of contrast media is required, withdraw the guide wire and inject the contrast media from the catheter hub with a syringe.

**CAUTION** Slowly inject a small amount of contrast media and observe the flow from the distal end of the product to ensure the lumen is maintained. When the product is occluded, it may result in deformation or damage of the product.

### 3. Removal of the product

3-1 Open the hemostasis valve.

3-2 Withdraw the product along the guide wire leaving the guide wire inside the vessel.

3-3 Close the hemostasis valve after removing the product.

#### CAUTIONS

- Confirm position of the guide wire under fluoroscopy when the product is removed.
- Remove the product, the guide wire and the guiding catheter altogether if any resistance is felt while withdrawing the product.
- Rinse residual blood on the surface of the product in a heparinized saline solution bath. When it is difficult to remove the residual blood, carefully wipe the surface with gauze moistened with heparinized saline solution. Flush the inner lumen of the product in order to remove residual blood inside the lumen.
- Carefully handle the product not to kink it. If it gets kinked, stop using it. Using the kinked product may result in separation or breakage of the product.
- Take care when re-inserting an angled guide wire into the product which remains in the coronary artery. Just before the tip of the guide wire protrudes from the tip of the product, stop advancing the guide wire and carefully withdraw the product so that the tip of the guide wire gradually protrudes from the product under high resolution fluoroscopy. The tip of the guide wire may spring out of the tip of the product, which may result in damage to blood vessels.

#### PRECAUTION FOR STORAGE

Avoid exposure to water, direct sunlight, extreme temperature, or high humidity during storage.

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© Terumo Corporation 2020-09
All brand names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners.