

# Glidewire Advantage Track™

Peripheral Guidewire

## LEVERAGE SUPERIOR PERFORMANCE

For your demanding CLI procedures

Available in 0.018" & 0.014"  
diameter sizes



The addition of **GLIDEWIRE ADVANTAGE TRACK™** Guidewire to your PAD/CLI toolkit provides you with **superior performance<sup>1,a</sup>** for your most challenging cases:

- **Superior Pushability:** High deformation resistance provides superior pushability to access and cross challenging above-the-knee (ATK) & below-the-knee (BTK) anatomy and lesions
- **Superior Trackability:** Provides superior ability to track and navigate anatomy of the vasculature
- **Best-in-Class Durability:** Provides superior shape retention for navigating small, tortuous vessels
- **Superior Lubricity:** TERUMO Glide Technology™ hydrophilic coating<sup>b</sup> provides superior lubricity to support ease of progression through vessels
- **Enhanced Visibility & Selectivity:** Gold coil tips enhance visibility in conjunction with a 25 cm distal taper designed to minimize vessel trauma and improve selectivity of vessels

<sup>a</sup> When compared to Abbott Hi-Torque Command 18 LT and Boston Scientific V-18

<sup>b</sup> Referring to the distal portion of GLIDEWIRE ADVANTAGE TRACK™ wire which is GLIDEWIRE®

# PUSHING BOUNDARIES

Terumo Interventional Systems is **committed to your success** with innovative procedural solutions and ongoing support for your most challenging cases.


We are relentlessly seeking new ways to help you achieve **better outcomes for more patients.**



## GLIDEWIRE ADVANTAGE TRACK™ Peripheral Guidewire

PRODUCT CODE	DIAMETER	TOTAL LENGTH	SPIRAL PTFE PROXIMAL COATING LENGTH	DISTAL GLIDEWIRE® LENGTH	FLEXIBLE TAPER LENGTH (TIP)	TIP SHAPE
GAT1418	0.014"	180 cm	155 cm	25 cm	1 cm	35° Angle
GAT1430	0.014"	300 cm	275 cm	25 cm	1 cm	35° Angle
GAT1818	0.018"	180 cm	155 cm	25 cm	1 cm	35° Angle
GAT1830	0.018"	300 cm	275 cm	25 cm	1 cm	35° Angle

### Technical Specifications



LENGTHS	DIAMETERS	FLEXIBLE TAPER LENGTH	TIP SHAPE	DISTAL MARKER	OTHER
180 cm 300 cm	0.014" 0.018"	1 cm	35° Angle	2 cm gold coil	<ul style="list-style-type: none"> <li>25 cm GLIDEWIRE® distal end</li> <li>Spiral PTFE-coated stainless steel core</li> </ul>

**FIND OUT MORE** ☎ Phone: 800.862.4143 🌐 [terumo.com](http://terumo.com) 📠 Fax: 800.411.5870

#### Reference:

1. Data on file.

**RX ONLY.** Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

©2020 Terumo Medical Corporation. All rights reserved. The third-party trademarks used herein are trademarks of their respective owners. All other brand names are trademarks or registered trademarks of Terumo. PM-03199

