

PINNACLE Introducer Sheath (Standard)			
SIZE (Fr)	SHEATH OD (mm / in)	TIP ID (mm / in)	SHEATH ID (mm / in)
4	1.96 / 0.077	1.40 / 0.055	1.56 / 0.061
5	2.29 / 0.090	1.73 / 0.068	1.89 / 0.074
6	2.62 / 0.103	2.06 / 0.081	2.22 / 0.087
7	2.95 / 0.116	2.40 / 0.094	2.55 / 0.100
8	3.32 / 0.131	2.73 / 0.107	2.92 / 0.115
9	3.73 / 0.147	3.09 / 0.122	3.33 / 0.131
10	4.39 / 0.173	3.46 / 0.136	3.69 / 0.145
11	4.83 / 0.190	3.82 / 0.150	4.13 / 0.163

0.035" for 4 Fr; 0.038" mini stainless steel spring coil guidewire.

PINNACLE PRECISION Introducer Sheath			
SIZE (Fr)	SHEATH OD (mm / in)	TIP ID (mm / in)	SHEATH ID (mm / in)
4	1.96 / 0.077	1.40 / 0.055	1.56 / 0.061
5	2.29 / 0.090	1.73 / 0.068	1.89 / 0.074
6	2.62 / 0.103	2.06 / 0.081	2.22 / 0.087
7	2.95 / 0.116	2.40 / 0.094	2.55 / 0.100
8	3.32 / 0.131	2.73 / 0.107	2.92 / 0.115

0.021" mini nitinol or stainless steel guidewire 21G to 19G tapered needle.

PINNACLE TIF TIP Introducer Sheath			
SIZE (Fr)	SHEATH OD (mm / in)	TIP ID (mm / in)	SHEATH ID (mm / in)
5	2.29 / 0.090	1.73 / 0.068	1.89 / 0.074
6	2.62 / 0.103	2.06 / 0.081	2.22 / 0.087
7	2.95 / 0.116	2.40 / 0.094	2.55 / 0.100
8	3.32 / 0.131	2.73 / 0.107	2.92 / 0.115

0.035" or 0.038" mini stainless steel spring coil guidewire.

PINNACLE R/O II Radiopaque Marker Introducer Sheath			
SIZE (Fr)	SHEATH OD (mm / in)	TIP ID (mm / in)	SHEATH ID (mm / in)
4	2.11 / 0.083	1.40 / 0.055	1.55 / 0.061
5	2.44 / 0.096	1.73 / 0.068	1.88 / 0.074
6	2.77 / 0.109	2.06 / 0.081	2.21 / 0.087
7	3.10 / 0.122	2.40 / 0.094	2.54 / 0.100
8	3.48 / 0.137	2.73 / 0.107	2.92 / 0.115
9	3.89 / 0.153	3.09 / 0.122	3.33 / 0.131

0.035" for 4 Fr; 0.038" mini stainless steel spring coil guidewire.

PINNACLE R/O II HiFlo Introducer Sheath			
SIZE (Fr)	SHEATH OD (mm / in)	TIP ID (mm / in)	SHEATH ID (mm / in)
6	2.77 / 0.109	2.06 / 0.081	2.21 / 0.087
7	3.10 / 0.122	2.40 / 0.094	2.54 / 0.100
8	3.48 / 0.137	2.73 / 0.107	2.92 / 0.115

0.035" Optional mini stainless steel guidewire.

PINNACLE R/O II HiFlo Flow Rates		
SIZE (Fr)	WITH 3-Way Stopcock	WITHOUT 3-Way Stopcock
6	377 ml / min	486 ml / min
7	416 ml / min	610 ml / min
8	478 ml / min	698 ml / min

Sheath French Size Color Indicator

4Fr	5Fr	6Fr	7Fr
8Fr	9Fr	10Fr	11Fr

FIND OUT MORE

☎ Phone: 800.862.4143 🌐 terumo.com 📠 Fax: 800.411.5870

RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.

Source: Data on File.

©2019 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of Terumo. PM-02563