## Glidecath<sup>®</sup>

**Hydrophilic Coated Catheter** 



## **Dead Space Volume Measurements**<sup>1</sup>

	Fr	L [mm]	V [mm3]	V [mL]
GLIDECATH® 4 & 5 Fr	4	600	519.54	0.52
	4	650	562.84	0.56
	4	700	606.13	0.61
	4	800	692.72	0.69
	4	1000	865.90	0.87
	4	1200	1039.08	1.04
	4	1250	1082.38	1.08
	4	1350	1168.97	1.17
	4	1500	1298.85	1.30
	5	650	640.38	0.64
	5	700	689.64	0.69
	5	800	788.16	0.79
	5	1000	985.20	0.99
	5	1100	1083.72	1.08
GLIDECATH® 5 Fr XP	5	650	759.84	0.76
	5	800	935.19	0.94
	5	1000	1168.99	1.17
	5	1250	1461.23	1.46
	5	1300	1519.68	1.52
	5	1350	1578.13	1.58
	5	1500	1753.48	1.75

FIND OUT MORE US: 800.862.4143 CA: 833.883.7866



**US:** terumois.com CA: terumocanada.ca

Reference: 1. Data on file.

RX ONLY. The advertisement is directed to physicians only, and not to consumers. Refer to product labels and packaging insert for complete warnings, precautions, potential complications, and instructions for use. Products may not have regulatory approval in all countries. Please contact your local sales representative if you have questions about the availability of products in your area.

